



iCare HOME 1



iCare HOME 2

FEATURES COMPARISON

2 CR123A batteries	Charging Device	4 AA batteries
Compact rectangle	Device Design	Contoured and compact, elegant modern design.
Insert probe directly from container. Probe is ejected when device turns off.	Probe Loading & Stability	Probe applicator insertion.. Probe stays in device even when off
Easy to operate.	Forehead And Cheek Rest	Adjust directly at support vs dial, stays securely in place, even in carrying case.
Symmetrical green ring, audio feedback (beeps)	Determining Forehead & Cheek Rest Distance	New smart light guide offers further improved usability
Blinking red light, audio cues, icon lit on display panel	Error Codes	Guidance on interactive display panel, blinking red light, audio cues
Wired connection to Android with iCare Patient app or PC with EXPORT	Measurements Transfer	Bluetooth or wired connection to iCare Patient 2 Apple/Android app or PC with EXPORT
View in iCare CLINIC web portal after measurements have been transferred	Measurements Viewing	View in iCare Patient 2 Apple/Android app or in iCare CLINIC web portal after measurements have been transferred
Upright	Body Position	Supine, reclined, or upright position

PRICING

PURCHASE		PURCHASE	
MSRP	OUR PRICE	MSRP	INTRO PRICE
\$2,995	\$2,399	\$3,195	\$2,945
RENTAL		RENTAL	
1 st WEEK	ADDITIONAL WEEKS	1 st WEEK	ADDITIONAL WEEKS
\$250	\$175	\$275	\$195

ENLIVENED INCLUDES FREE:

- 1-on-1 training and educational video conference session
 - Getting started instructions
 - Unlimited email support



We're Always Here for You

Just contact a member of our Support Team by email at support@enlivened.com or call **1-844-844-4200**

We Empower **You** to Create Health